

Product Insert

CYCLOSERINE-CEFOXITIN FRUCTOSE AGAR (CCFA)

Products

AS-213 Cycloserine-Cefoxitin Fructose Agar (CCFA) 1 plate / pkg

Intended Use

Cycloserine-Cefoxitin Fructose Agar (CCFA) is an enriched selective and differential media used for the isolation and presumptive identification of *Clostridium difficile*, a recognized cause of pseudomembranous (antimicrobial agent-associated) colitis.

Summary

CCFA is an enriched selective and differential media used for the isolation of *Clostridium difficile*. The basic nutritive base consists of animal peptones and fructose, and is supplemented with cefoxitin and cycloserine at concentrations that inhibit the growth of most normal fecal flora. Cycloserine will inhibit gram-negative bacteria, while cefoxitin will inhibit both gram-positive and gram-negative organisms. *Clostridium difficile* is not inhibited on CCFA, and when grown, will exhibit a characteristic yellow, ground-glass colonial morphology. Neutral red is added as a pH indicator. The presence of *C. difficile* will turn the indicator from pink/orange to yellow due to the amino acids being utilized by the organism, increasing the overall pH causing the color change. This media is prepared, stored, and dispensed under oxygen-free conditions to prevent the formation of oxidized products prior to use.

Formulation

Proteose Peptone	40.00	g
Sodium Phosphate Dibasic	5.00	g
Potassium Phosphate Monobasic	1.00	g
Sodium Chloride	2.00	g
Magnesium Sulfate Heptahydrate	0.20	g
Fructose	6.00	g
Agar	15.00	g
Neutral Red (0.3% solution)	3.00	mL
Cycloserine (10.0% solution)	5.00	mL
Cefoxitin (1.56% solution)	1.00	mL
DI Water	1.00	L

Final pH: 7.2 ± 0.2 at 25 °C

Final weight: 16.0 g ± 1.6 g

*Approximate formula. Adjusted and/or supplemented as required to meet performance criteria.

Precautions

For *IN VITRO DIAGNOSTIC USE* only. Utilize approved biohazard precautions and aseptic technique when using this product. This product is for use only by properly trained and qualified personnel. Sterilize all biohazard waste prior to disposal.

Storage and Shelf Life

Storage: Upon receipt, store at 2 – 8°C in original package until used. Avoid overheating or freezing. Do not use media if there are signs of deterioration (shrinking, cracking, or discoloration due to oxidation of media) or contamination. The expiration date applies to the product in its original packaging and stored as directed. Do not use product past the expiration date shown on the container.

Shelf Life: 90 days from date of manufacture.

Procedure

Specimen Collection: Protect specimens for anaerobic culture from oxygen during collection, transportation, and processing. Consult appropriate references for detailed instructions concerning collection and transportation of anaerobes.

Methods for Use: CCFA should be inoculated directly with clinical specimen or from a broth that has been previously inoculated from a clinical specimen. Streak plates with inoculum to obtain isolated colonies and immediately place in an anaerobic atmosphere, incubating at 35-37°C for 18-48 hours. Quantitation of *C. difficile* in a specimen may be clinically useful, which can be achieved by thoroughly mixing a serial 10-fold dilution of the specimen in an anaerobic environment followed by plating the dilutions onto CCFA media. Detailed instructions for processing anaerobic cultures can be found in the appropriate references.

Materials Required, But Not Provided

Standard microbiological supplies and equipment such as loops, saline blanks, slides, staining supplies, microscope, incinerator / autoclave, incubators, anaerobic chamber / anaerobic jars, other culture media, and serological / biochemical reagents.

Interpretation of Results

If used properly, CCFA supports good growth of *C. difficile*. After 24-48 hours, most colonies of *C. difficile* are large, circular, and yellow, with the yellow coloration extending 2–3 mm beyond the colony into the initially pink/orange media. These colonies fluoresce golden yellow/chartreuse under long-wavelength UV light. Most other bacteria are inhibited on this media. At 48 hours, colonies of most other organisms (e.g. *Lactobacilli* and yeast), which may grow, are very small (pinpoint to 0.5 mm in diameter) and do not fluoresce golden-yellow.

Limitations

CCFA will not provide complete information for identification of bacterial isolates. Rare strains of *C. difficile* may be inhibited. Plates must be examined no later than after 48 hours of incubation for optimal selectivity, as after 3 – 5 days of incubation, significant numbers of colonies other than *C. difficile* may grow. A test for aerotolerance should be used to confirm that each colony type is an obligate anaerobe. Consult reference materials for additional information.

Quality Control

The following organisms are routinely used for quality control performance testing at Anaerobe Systems.

Organism Tested	ATCC #	Results	Time	Special Reaction
<i>Bacteroides fragilis</i>	25285	No Growth		
<i>Prevotella melaninogenica</i>	25845	No Growth		
<i>Fusobacterium necrophorum</i>	25286	No Growth		
<i>Fusobacterium nucleatum</i>	25586	No Growth		
<i>Clostridium perfringens</i>	13124	No Growth		
<i>Peptostreptococcus anaerobius</i>	27337	No Growth		
<i>Staphylococcus aureus</i> or <i>Enterococcus faecalis</i>	25923 29212	No Growth		
<i>Escherichia coli</i>	25922	No Growth		
<i>Clostridium difficile</i>	9689	Growth	24 hrs	Yellow coloration, chartreuse fluorescence

User Quality Control: The final determination to the extent and quantity of user laboratory quality control must be determined by the end user.

If sterility testing is to be performed on this product, a representative sample of the lot(s) should be incubated anaerobically and aerobically for 48 – 96 hours.

If the nutritive/inhibitory capacity of this media is to be tested for performance, it is recommended that the following ATCC organisms be evaluated for growth/inhibition.

Organism	ATCC #	Expected Growth	Special Reaction
B. fragilis	25285	No Growth	
E. coli	25922	No Growth	
S. aureus	25923	No Growth	
C. difficile	9689	24 hours	Yellow coloration, chartreuse fluorescence

Physical Appearance: CCFA should appear translucent pink/orange in color.

References

1. Dowell, V. R., Jr. and T. M. Hawkins. 1987. *Laboratory Methods in Anaerobic Bacteriology*. CDC Laboratory Manual. USDHHS CDC. Atlanta, GA 30333.
2. Dowell, V. R., Jr. and G. L. Lombard. 1981. *Presumptive Identification of Anaerobic Non-sporeforming Gram-negative Bacilli*. USDHEW, CDC. Atlanta, GA 30333.
3. Dowell, V. R., Jr., G. L. Lombard, F. S. Thompson and A. Y. Armfield. 1977. *Media for the Isolation, Characterization, and Identification of Obligately Anaerobic Bacteria*. USDHHS, CDC, Atlanta, GA 30333.
4. Holdeman, L. V., F. P. Cato and W. E. C. Moore. 1987. *Anaerobe Laboratory Manual*. Virginia Polytechnic Institute and State University. Blacksburg, VA 24061.
5. Jousimies-Somer, H. R., P. Summanen, D. M. Citron, E. J. Baron, H. M. Wexler and S. M. Finegold. 2002. *Wadsworth – KTL Anaerobic Bacteriology Manual*. Star Publishing Co., Belmont, CA 94002.
6. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard- Third Edition*. (2004). CLSI document M22-A3. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898.
7. George, W. L., V.L. Sutter, D. Citron, S. Finegold. 1979. Selective and Differential Medium for Isolation of *Clostridium difficile*. *Journal of Clinical Microbiology*. 9:214-219.

Revision Date: 09/27/17